Biacore™ 1 シリーズ

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Systems specifications

Technical specifications and characteristics

| | |
|----------------------------|--|
| Detection technology | Surface plasmon resonance (SPR) biosensor |
| Information provided | Kinetic and affinity data (k_a , k_d , K_D), specificity, selectivity, screening data, epitope binning, concentration and relative potency data |
| Data presentation | Monitoring of real-time sensorgrams or evaluation data for result tables and result plots |
| Analysis time per cycle | Typically 2 to 15 min |
| Automation | 60 h unattended run time for Biacore™ 1K 72 h unattended run time for Biacore™ 1K+ and Biacore™ 1S+ |
| Sample type | Small molecule drug candidates to high molecular weight proteins (also DNA, RNA, polysaccharides, lipids, cells, and viruses) in various sample environments (e.g., in DMSO-containing buffers, plasma, and serum) |
| Required sample volume | Injection volume plus 20 to 40 μ L (application-dependent) |
| Injection volume | 1 to 400 μ L |
| Flow rate range | 1 to 100 μ L/min |
| Flow cell volume | 60 nL |
| Flow cell height | 50 μ m |
| Data collection rate | 1 or 10 Hz for Biacore™ 1K and Biacore™ 1K+ 1, 10 Hz or 40 Hz for Biacore™ 1S+ |
| Sample/reagent capacity | 1 \times 96- or 384-well microplate, normal, and deep-well 1 \times reagent rack with 21-43 positions compatible with 0.7-4.4 mL vials (Biacore™ 1K) 2 \times 96- or 384-well microplates, normal, and deep-well 2 \times reagent racks with 21-43 positions compatible with 0.7-4.4 mL vials (Biacore™ 1K+ and Biacore™ 1S+) |
| Typical run times | Clean screen (384-well plate): 6 h Binding level screen (384-well plate): 15 h Affinity screen (384- samples): 27 h Kinetic analysis (30 samples): 15 h Concentration analysis (24 samples): 3h Epitope binning, 5 \times 5 array (25 samples): 1.5 h |
| Analysis temperature range | 25°C to 37°C for Biacore™ 1K and Biacore™ 1K+ 4°C to 40°C for Biacore™ 1S+ (at least 20°C below ambient temperature) |
| Sample storage | 4°C to 37°C for Biacore™ 1K and Biacore™ 1K+ (at least 18°C below ambient temperature) 4°C to 40°C for Biacore™ 1S+ (at least 18°C below ambient temperature) |

| | |
|--------------------------------------|--|
| Sample refractive index range | 1.33 to 1.39 |
| In-line reference subtraction | Automatic |
| Number of flow cells | 6 in 1 channel |
| Dimensions (W \times H \times D) | 755 \times 725 \times 666 mm |
| Net weight total | 95 kg (Biacore™ 1K) 96 kg (Biacore™ 1K+ and Biacore™ 1S+) |
| Mains requirements | Processing unit: Autorange voltage 100 to 240 V~, frequency 50/60 Hz |
| Power consumption | Processing unit: max. 350 VA |

Minimum computer requirements

| |
|---|
| CPU with at least four cores |
| At least 16 GB internal memory |
| At least 200 GB free hard disk space |
| Screen resolution at least 1920 \times 1080 |
| One USB2 port available for instrument connection |

Minimum network SQL Server requirements

| |
|---|
| See recommendations from Microsoft® regarding hardware and operating systems. |
| Uninterrupted power supply |

Typical working ranges

| | |
|--------------------------------------|---|
| Association rate constant (k_a) | Proteins: up to $3 \times 10^9 \text{ M}^{-1} \text{ s}^{-1}$ LMW molecules: up to $5 \times 10^7 \text{ M}^{-1} \text{ s}^{-1}$ |
| Dissociation rate constant (k_d) | 10^{-6} to 1 s^{-1} for Biacore™ 1K and Biacore™ 1K+ 10^{-6} to 6 s^{-1} for Biacore™ 1S+ |
| Sample concentration | ≥ 1 picomolar (pM) |
| Molecular weight detection | No lower limit for organic molecules |
| Short term noise typically | < 0.03 RU (RMS) for Biacore™ 1K and Biacore™ 1K+ < 0.01 RU (RMS) for Biacore™ 1S+ |
| Baseline drift typically | < 0.3 RU/min |
| Blank subtracted drift | < +/-0.003 RU/min |
| Immobilized interactant consumption | Typically 0.03 to 3 μ g/flow cell |

Data handling and storage

| | |
|---------------------|---|
| Operating system | Windows® 10 (Professional or Enterprise), 64-bit, English version Note: The functionality of Biacore™ Insight Software and Biacore™ system is verified using an English version of Windows®. Other languages can cause issues. |
| Interfacing | Import of sample data and export of results possible |
| Licenses | Multiple licenses available |
| Server requirements | Includes SQL Server Express 2019. Security and performance improvements are seen with SQL Server Standard or SQL Server Enterprise 2017, or 2019 (available separately from Microsoft®) |

Notes: The server is supplied by the end user. Contact your local representative for the latest information regarding on-site requirements.

Compliance

| | |
|-------------------------------------|--|
| Compliant with | CE, cETLus, EAC, FCC, ICES-001, KC, RCM, UKCA |
| Safety | IEC/EN/UL/CSA-C22.2 61010-1, IEC/EN/UL/CSA-C22.2 61010-2-081, EN ISO 12100 |
| Electromagnetic compatibility (EMC) | EN/IEC 61326-1, FCC Part 15 B, ICES-001 |
| Environmental | EN 63000, China RoHS |

For local office contact information, visit cytiva.com/contact

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. Biacore, Biacore Multi Cycle Kinetics (MCK), and Biacore Single Cycle Kinetics (SCK) are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Excel, Microsoft, PowerPoint, and Windows are registered trademarks of the Microsoft Corporation. Any other third party trademarks are the property of their respective owners.

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© 2022 Cytiva