Biacore™ 1 シリーズ

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Systems specifications

Technical specifications and characteristics

Detection technology	Surface plasmon resonance (SPR) biosensor
Information provided	Kinetic and affinity data (k_a , k_d , K_D), specificity, selectivity, screening data, epitope binning, concentration and relative potency data
Data presentation	Monitoring of real-time sensorgrams or evaluation data for result tables and result plots
Analysis time per cycle	Typically 2 to 15 min
Automation	60 h unattended run time for Biacore™ 1K 72 h unattended run time for Biacore™ 1K+ and Biacore™ 1S+
Sample type	Small molecule drug candidates to high molecular weight proteins (also DNA, RNA, polysaccharides, lipids, cells, and viruses) in various sample environments (e.g., in DMSO-containing buffers, plasma, and serum)
Required sample volume	Injection volume plus 20 to 40 μ L (application-dependent)
Injection volume	1 to 400 μ L
Flow rate range	1 to 100 μ L/min
Flow cell volume	60 nL
Flow cell height	50 μ m
Data collection rate	1 or 10 Hz for Biacore™ 1K and Biacore™ 1K+ 1, 10 Hz or 40 Hz for Biacore™ 1S+
Sample/reagent capacity	1 \times 96- or 384-well microplate, normal, and deep-well 1 \times reagent rack with 21-43 positions compatible with 0.7-4.4 mL vials (Biacore™ 1K) 2 \times 96- or 384-well microplates, normal, and deep-well 2 \times reagent racks with 21-43 positions compatible with 0.7-4.4 mL vials (Biacore™ 1K+ and Biacore™ 1S+)
Typical run times	Clean screen (384-well plate): 6 h Binding level screen (384-well plate): 15 h Affinity screen (384- samples): 27 h Kinetic analysis (30 samples): 15 h Concentration analysis (24 samples): 3h Epitope binning, 5 \times 5 array (25 samples): 1.5 h
Analysis temperature range	25°C to 37°C for Biacore™ 1K and Biacore™ 1K+ 4°C to 40°C for Biacore™ 1S+ (at least 20°C below ambient temperature)
Sample storage	4°C to 37°C for Biacore™ 1K and Biacore™ 1K+ (at least 18°C below ambient temperature) 4°C to 40°C for Biacore™ 1S+ (at least 18°C below ambient temperature)

Sample refractive index range	1.33 to 1.39
In-line reference subtraction	Automatic
Number of flow cells	6 in 1 channel
Dimensions (W \times H \times D)	755 \times 725 \times 666 mm
Net weight total	95 kg (Biacore™ 1K) 96 kg (Biacore™ 1K+ and Biacore™ 1S+)
Mains requirements	Processing unit: Autorange voltage 100 to 240 V~, frequency 50/60 Hz
Power consumption	Processing unit: max. 350 VA

Minimum computer requirements

CPU with at least four cores
At least 16 GB internal memory
At least 200 GB free hard disk space
Screen resolution at least 1920 \times 1080
One USB2 port available for instrument connection

Minimum network SQL Server requirements

See recommendations from Microsoft® regarding hardware and operating systems.
Uninterrupted power supply

Typical working ranges

Association rate constant (k_a)	Proteins: up to $3 \times 10^9 \text{ M}^{-1} \text{ s}^{-1}$ LMW molecules: up to $5 \times 10^7 \text{ M}^{-1} \text{ s}^{-1}$
Dissociation rate constant (k_d)	10^{-6} to 1 s^{-1} for Biacore™ 1K and Biacore™ 1K+ 10^{-6} to 6 s^{-1} for Biacore™ 1S+
Sample concentration	≥ 1 picomolar (pM)
Molecular weight detection	No lower limit for organic molecules
Short term noise typically	< 0.03 RU (RMS) for Biacore™ 1K and Biacore™ 1K+ < 0.01 RU (RMS) for Biacore™ 1S+
Baseline drift typically	< 0.3 RU/min
Blank subtracted drift	< +/-0.003 RU/min
Immobilized interactant consumption	Typically 0.03 to 3 μ g/flow cell

Data handling and storage

Operating system	Windows® 10 (Professional or Enterprise), 64-bit, English version Note: The functionality of Biacore™ Insight Software and Biacore™ system is verified using an English version of Windows®. Other languages can cause issues.
Interfacing	Import of sample data and export of results possible
Licenses	Multiple licenses available
Server requirements	Includes SQL Server Express 2019. Security and performance improvements are seen with SQL Server Standard or SQL Server Enterprise 2017, or 2019 (available separately from Microsoft®)

Notes: The server is supplied by the end user. Contact your local representative for the latest information regarding on-site requirements.

Compliance

Compliant with	CE, cETLus, EAC, FCC, ICES-001, KC, RCM, UKCA
Safety	IEC/EN/UL/CSA-C22.2 61010-1, IEC/EN/UL/CSA-C22.2 61010-2-081, EN ISO 12100
Electromagnetic compatibility (EMC)	EN/IEC 61326-1, FCC Part 15 B, ICES-001
Environmental	EN 63000, China RoHS

Ordering information

Product	Product code
Biacore™ 1K system Includes: Biacore™ 1K instrument (29712576); 2 licenses for Biacore™ Insight Control and Evaluation Software (29310602)	29726017
Biacore™ 1K+ system Includes: Biacore™ 1K+ instrument (29712577); 2 licenses for Biacore™ Insight Control and Evaluation Software (29310602)	29726018
Biacore™ 1S+ system Includes: Biacore™ 1S+ instrument (29712578); 2 licenses for Biacore™ Insight Control and Evaluation Software (29310602)	29726019
Biacore™ Insight Software extensions	Various licenses ¹
Biacore™ 1K+ Upgrade Kit Includes: Upgrade of Biacore™ 1K to Biacore™ 1K+ performed by Cytiva field service	29719933

¹ See cytiva.com/biacore for details of the various e-licenses available.

For local office contact information, visit cytiva.com/contact

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. Biacore, Biacore Multi Cycle Kinetics (MCK), and Biacore Single Cycle Kinetics (SCK) are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Excel, Microsoft, PowerPoint, and Windows are registered trademarks of the Microsoft Corporation. Any other third party trademarks are the property of their respective owners.

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© 2022 Cytiva