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
Show by technique **Affinity**

Column type **Any**

Show only suggested columns Column Properties...

Column volume **0.400** ml

Pressure limit pre-column **1.00** MPa [0.02 - 20.00]

Pressure limit delta-column **1.00** MPa [0.02 - 20.00]

Use flow restrictor

Column position **By-pass**

Flow rate **16.000** ml/min [0.000 - 25.000]

Control the flow to avoid overpressure

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0.00 Noise reduction UV: (0.2)#UV averaging time (sec)

Equilibration

Reset UV monitor (recommended if the equilibration occurs before the purification).

Use the same flow rate as in Method Settings Use the same inlets as in Method Settings

Flow rate **16.000** ml/min [0.000 - 25.000]

Inlet A **A1**

Inlet B **B1** **0.0** % B [0.0 - 100.0]

Fill the system with the selected buffer

Equilibrate until

the total volume is **6.00** ml

the following condition is met

Conductivity greater than

Conductivity greater than **0.00** mS/cm [0.00 - 1000.00]

Accepted pH fluctuation **0.10** [0.00 - 14.00]

Accepted UV fluctuation **0.10** mAU [0.00 - 6000.00]

Accepted conductivity fluctuation **0.10** mS/cm [0.00 - 300.00]

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Sample Application

Use the same flow rate as in Method Settings
Flow rate ml/min [0.000 - 50.000]

Inject sample from loop
 Inject sample directly onto column

Sample inlet

Inject fixed sample volume ml
 Inject all sample using air sensor
 Set maximum volume to ml
 Finalize sample injection ml

Wash sample pump with buffer
 Prime sample inlet with ml
 Wash sample pump with buffer after sample application

Interrupt sample application at UV mAU [-6000.0 - 6000.0]

Fractionate

in waste (do not collect)
 using outlet valve
 using fraction collector
Fraction collector 1

Fractionation settings

Fractionation type
Fractionation destination
Peak fractionation destination
Fixed fractionation volume ml [0.01 - 20000.00]
Peak fractionation volume ml [0.01 - 20000.00]

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Fibro Unit Wash

Use the same flow rate as in Method Settings
Flow rate ml/min [0.000 - 25.000]

Use the same inlets as in Method Settings
Inlet A
Inlet B % B [0.0 - 100.0]
 Fill the system with the selected buffer

Wash until

the total volume is ml
 the following condition is met

Stable UV mAU [-6000.0 - 6000.0]
UV less than mAU [-6000.0 - 6000.0]
Stability time min [0.02 - 1000.00]
Accepted UV fluctuation mAU [0.00 - 6000.00]
Maximum wash volume ml

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Elution

Use the same flow rate as in Method Settings Use the same inlets as in Method Settings
 Flow rate: ml/min [0.000 - 25.000] Inlet A: Inlet B:

Up flow

Isocratic elution
 Volume: ml % B [0.0 - 100.0] Fill the system with the selected buffer

Gradient elution
 Start at: % B [0.0 - 100.0] Fill the system with the selected buffer

Type	Target %B (0-100)	Length (ml)
1 Step with fill	100.0	2.00

Note: A gradient delay is automatically added, provided that the last gradient segment is linear

Fractionate: in waste (do not collect) using outlet valve using fraction collector

Fractionation settings:
 Fractionation type:
 Fractionation destination:
 Peak fractionation destination:
 Fixed fractionation volume: ml [0.01 - 20000.00]
 Peak fractionation volume: ml [0.01 - 20000.00]

Start fractionation after ml (only for isocratic elution)

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Re Equilibration	
<input type="checkbox"/> Reset UV monitor (recommended if the equilibration occurs before the purification).	
<input checked="" type="checkbox"/> Use the same flow rate as in Method Settings Flow rate <input type="text" value="16.000"/> ml/min [0.000 - 25.000]	<input checked="" type="checkbox"/> Use the same inlets as in Method Settings Inlet A <input type="text" value="A1"/> Inlet B <input type="text" value="B1"/> <input type="text" value="0.0"/> % B [0.0 - 100.0]
<input type="checkbox"/> <u>Fill the system with the selected buffer</u>	
Equilibrate until <input checked="" type="radio"/> the total volume is <input type="text" value="10.00"/> ml <input type="radio"/> the following condition is met	
Conductivity greater than <input type="text" value="0.00"/> mS/cm [0.00 - 1000.00]	
Accepted pH fluctuation <input type="text" value="0.10"/> [0.00 - 14.00]	
Accepted UV fluctuation <input type="text" value="0.10"/> mAU [0.00 - 6000.00]	
Accepted conductivity fluctuation <input type="text" value="0.10"/> mS/cm [0.00 - 300.00]	
Signal stable for <input type="text" value="1.00"/> min [0.02 - 1000.00]	
Maximum equilibration volume <input type="text" value="4.00"/> ml	

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