
Cytiva Webinar

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
Sensor chip	抗体種
Protein A	Human IgG1, IgG2 and IgG4 mammalian IgG Stability of mouse IgG1 is quite low, so kinetic studies challenging
Protein G	Human, rat, mouse, rabbit, sheep, goat, pig, cow, horse IgG Not bind chicken IgG, human IgA, IgD, IgE and IgM
Protein L (for yes/no screening)	Kappa light chain subtypes 1,3 and 4 Not bind bovine immunoglobulins (so useful for antibodies from serum based culture)

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



Assay setup

	Report point-based screening	Sensorgram-based screening
特徴	Blank が不要 <ul style="list-style-type: none"> • Capture level • Stability early, Stability late 	要Blank <ul style="list-style-type: none"> • Off-rate ranking • Kinetics screening
(Analyte) Flow rate	> 30 uL/min	> 30 uL/min
(Analyte) Injection type	Low sample consumption	High performance
(Analyte) Contact time	60 sec	60-120 sec
(Analyte) Dissociation time	120 sec	180 sec
(Analyte) Pooling	Optional	Optional
(Analyte) Concentration	Use the same conc. In all cycles	Use the same conc. In all cycles
(Ligand) Flow rate	10 uL/min	10 uL/min
(Ligand) Contact time	60-180 sec	60-180 sec

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	Ligand	Analyte
Cycle 04	●	○ Blank
Cycle 05	●	●
Cycle 06	●	○
Cycle 07	●	●
Cycle 08	●	○
Cycle 09	●	●
Cycle 10	●	○
Cycle 11	●	●
Cycle 12	●	○
Cycle 13	●	●
Cycle 14	●	○
Cycle 15	●	●

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	Ligand	Analyte
Cycle 04	● Rep	○ Blank
Cycle 05	●	●
Cycle 06	●	●
Cycle 07	●	●
Cycle 08	●	●
Cycle 09	●	●
Cycle 10	● Rep	○ Blank
Cycle 11	●	●
Cycle 12	●	●
Cycle 13	●	●
Cycle 14	●	●
Cycle 15	●	●

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Early kinetic screening of hybridomas for confident antibody selection using Biacore A100

- Rapid, kinetic screening of mAbs direct from hybridoma samples
 - efficient screening rate of ~800 hybridomas/day
 - effective selection of mAbs with desired kinetic profiles
- High resolution kinetic characterization of selected mAbs
 - affinity and kinetic parameters determined directly from crude samples
 - dedicated evaluation software processes kinetic data in minutes
 - selection of candidates based on therapeutically relevant binding properties
 - parallel analysis under different buffer conditions for characterization of stability of antibody performance
- More confident mAb selection earlier in the development process compared to end-point assays such as ELISA


Figure 1. Biacore A100 offers rapid kinetic screening and high resolution kinetic characterization of mAbs direct from crude hybridoma samples in

Evaluation of screening results using report point values

Label-free interaction analysis produces a real-time plot of binding response against time (the *sensorgram*) over the entire course of binding and dissociation. Binding responses at specific selected times (report points) can be taken, to rapidly plot binding characteristics of interest (see *Monitoring Interactions* at the end of this note). Figure 3 shows a plot of responses from the late association phase (*binding_late* report point) from each mAb-antigen interaction recorded in the screening assays. Since the *binding_late* response is related to the association rate of the interaction and *stability_late* is inversely related to the dissociation rate, these two values provide a good overview of the kinetic profile of each mAb.


Figure 3. Rapid kinetic screening of hybridomas, evaluated by report point plot (responses normalized to capture level). Areas corresponding to low/non-binders and highest affinity, as well as the 16 mAbs selected for further characterization are indicated on the figure.

As seen in Figure 3, mAbs from the 384 hybridoma samples displayed a wide range of kinetic profiles. A number of negative hybridomas were supplied blind by Abgenix, and all of these were determined to be non-binders by the kinetic screening assay (i.e. no known false-positives were detected). Using this plot, 16 hybridomas were selected for comprehensive kinetic characterization, covering a wide range of association and dissociation properties among the positive samples.

Estimation of affinity & kinetic parameters from rapid kinetic screen

Single binding curves can be used to derive approximate k_a and k_d values and hence, estimate affinities, although fully quantitative, high-resolution kinetic characterization requires analysis of multiple sensorgrams derived from a concentration series of analyte. To obtain an estimate of the affinity/kinetic parameter ranges corresponding to the profiles observed in the report point plots, sensorgrams were therefore fitted to a 1:1 interaction model and the approximate rate constants used to generate an off/on-rate map (Figure 4). From this log-scale plot of estimated k_d against estimated k_a , a generally good correlation in kinetic profiles was observed compared to those indicated by the report point plot (Figure 3). The off/on-rate map also showed that many mAbs lay within an approximate affinity range of 1-10 nM and that the three highest affinity examples lay in the sub-nanomolar range.


Figure 4. Rapid kinetic screening of hybridomas, evaluated by estimated rate constants. Log-scale plot of estimated k_d against estimated k_a . The isomeric lines corresponding to 10, 1, 0.1, 0.01 nM affinities are indicated by the red diagonals. Negative hybridomas and lower-level binders were excluded from this plot.

ent an important class whole biomedical research tools in the (as precision diagnostic reapeutics against a wide rances now enable the bodies directed against a) increasing demand for high quality binding data action.

lent, end-point Biacore™ systems for analysis of antibody based on functionally

App Note 79

Development of immunotherapeutics and immunization regimes

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Application Note 79 Biacore systems

Development of immunotherapeutics and immunization regimes

Characterization of serum antibody responses during development of an anti-IgE based treatment for allergy and asthma, using Biacore T100

- Rapid, high information content characterization of immune responses direct from rat and primate serum
 - serum response patterns monitored over time following multiple immunizations
 - IgG subclass rapidly determined
 - qualitative changes derived from kinetic profiles
- Serum response assays correlated well with observed biological effects
- Significant advantages shown over ELISA
 - higher quality data: much lower inter- and intra-assay variations

Biacore T100

Summary

- The recombinant IgE-derived immunotherapeutic proteins developed by Resistoria Pharmaceuticals AB show excellent immunogenic properties in the two animal systems tested, providing the basis for an effective immunotherapeutic approach to the treatment of allergy and asthma.
- Biacore T100 provided high information content characterization of anti-IgE responses over time and in relation to immunization dosing and adjuvance regimes
 - Increased anti-IgE responses correlated strongly with reductions in circulating IgE levels, demonstrating neutralizing activity of the immunotherapeutic protein-induced antibodies
 - Subclass analysis showed a clear correlation between subclass switching patterns and effects on circulating IgE levels among different individuals
 - Rapid kinetic profiling of primate anti-IgE responses provided information on the quality of the antibody response over time and in relation to the dosing regime.
- Response data showed a good general agreement with ELISA, but showed a much higher level of reproducibility.

Methods

IgE and anti-IgE ELISAs were carried out using standard methodologies at Resistoria Pharmaceuticals AB. Biacore T100 assays were performed using Series 5 Sensor Chip CM5 and an assay temperature of 37°C. Standard amine coupling procedures were used to immobilize a commercial rat IgE protein (rat studj) or a Cynomolgus monkey IgE-derived recombinant protein (primate studj) onto the sensor surface. An overview of the assay setups for the two animal studies is shown in Figure 8.

Figure 8. Overview of assay setups used in the immune response assays. A) Rat study: commercial rat IgE protein immobilized to detect anti-IgE Abs in immunized animals. Further serial injections of IgG subclass-specific Abs is identified subclass of antigen-bound serum anti-IgE Abs. B) Primate study: recombinant Cynomolgus monkey IgE immobilized on the sensor surface. C) Deviation of data from sensograms. Binding of proteins on the sensor surface generates an SPR signal that is monitored in real time over the whole binding interaction (total interaction time: 30 minutes in the example shown). Anti-IgE and subclass responses were determined from single report points, taken from specific time-points on the sensograms. Qualitative assessment of Abs was based on slopes taken from the early dissociation phase. In practice, the response level and dissociation profiling data were obtained from separate experiments, with sensor surface preparation and other conditions tailored for the specific assay.

Biacore

App Note 28-9777-72AA


Resolving a bottleneck in screening and characterization of recombinant antibody fragments using Biacore 4000

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Resolving a bottleneck in screening and characterization of recombinant antibody fragments using Biacore 4000

Antibodies and their derivatives are increasingly important as biotherapeutics, precision diagnostics, and essential tools for biological research. Techniques to screen and characterize large numbers of antibodies and provide confident selection of the best candidates for development are therefore essential. Biacore 4000 equipped with Biacore 4000 Antibody Extension Package has been used for efficient screening and characterization of antibody fragments for diagnostic purposes. Ten 96-well plates containing avian scFv fragments in crude cell lysates were screened and ranked in less than 20 h. A subset from this screen was selected for analysis with the 2-over-2 kinetics approach. This method

Results and discussion

Screening of scFv antibody fragments based on binding stability

A panel of 960 scFv fragments specific to human cardiac troponin I was screened using the assay setup described in Figure 2. Sensorgrams from 120 samples with the report points indicated as colored bars are shown in Figure 4.


Fig 4. Sensorgrams from 120 samples showing report points used for ranking.

Figure 5 shows a scatterplot of the report points stability_early plotted against stability_late. This type of plot can be used to rank binders based on binding stability. The figure shows two main groups of scFvs: non-binders close to the origin and binders higher along the diagonal. The best binders, with a high binding stability and slow dissociation, are shown as blue squares. In total, 960 samples were analyzed and ranked with respect to binding stability in 18 h. The data evaluation was accomplished in less than one hour.


Fig 5. Stability_early vs Stability_late plot for identification of stable binders (best binders shown as blue squares, green binders are not as stable).

Importantly, Biacore analysis of antibodies and their fragments allows both the expression levels of the clones and quality of the antibodies to be measured simultaneously. Figure 6 shows an example of how report point data from the scFv screen are plotted in Spotfire, providing an excellent overview of the performance of different clones.


Fig 6. Overview of the scFv screening. Antigen binding stability, stability_early (top) and stability_late (bottom), is plotted against capture level in spot 1 (left) and spot 3 (right). Wells without scFv (blue), wells with significant capture (red), and scFvs showing a stable antigen response (circles) are displayed. The green box indicates scFvs that are captured, but do not bind to the antigen. The best scFv clones show high capture level and high binding stability (blue circles).

The amount of captured scFv on spot 1 and 5 in the four flow cells is plotted on the x-axis, and indicates the capture level or expression level of the different clones. After scFv has been captured, antigen is injected over the spots. Response levels early and late in antigen dissociation are plotted on the y-axis to show the quality of the interaction of captured scFv with antigen. Uninteresting clones are distinguished from the good clones (close to the dashed lines), and the best clones show both high capture level and binding stability (slow off-rates). These are important criteria when selecting high quality antibodies and antibody fragments for further development, and were applied to the scFvs selected for 2-over-2 kinetic characterization.

Characterization with 2-over-2 kinetics

Obtaining kinetics data early in the development process enables better-informed selection since kinetic properties have consequences for both therapeutic and diagnostic agents. The 2-over-2 kinetics approach is particularly suitable for providing kinetic information already during screening, as the data can be obtained in half the time using significantly less reagents and sample volumes.

Evaluation was performed with the kinetic evolution tool in Biacore 4000 Evaluation Software, by combined kinetic fitting of the sensorgram data from spots 1-2 and 4-5, respectively. The evolution combines the different surface densities with the different analyte concentrations in one fit.

To investigate temperature stability of the binding, all samples have been analyzed at 25°C (Fig 7A) and 37°C (Fig 7B). Kinetic constants from the two temperatures were exported to Spotfire for further analysis. Temperature-stable, high-affinity scFvs were identified in Spotfire using on/off-rate maps (Fig 8).

1 x 10⁵ clones
 1 x 10² clones
 studied at each screening stage
 less for selecting recombinant antibodies
 for secondary screening, has the
 analysis bottleneck and provides high
 ic data.
 r (prostate-specific antigen, PSA)
 ostatic reagents, ideal antibodies
 and a high affinity that is

Assessing long-term glucose regulation by measurement of glycated hemoglobin using Biacore 8K

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Assessing long-term glucose regulation by measurement of glycated hemoglobin using Biacore™ 8K

Persistently elevated levels of blood glucose lead to irreversible attachment of glucose molecules to hemoglobin in red blood cells. The longer the hyperglycemia occurs, the more glucose binds to hemoglobin. By measuring the binding of antibodies that recognize glycated hemoglobin

Biacore systems has become an indispensable process standard in the early development of mAbs for diagnostic and pharmaceutical applications. In particular, kinetic antibody screening is established in the diagnostic antibody production process. Kinetic screening resolves

A Biacore Series S Sensor Chip CM5 was mounted into the Biacore 8K system and was normalized in instrument buffer (10 mM HEPES pH 7.4, 150 mM NaCl, 3 mM EDTA, 0.005% w/v Tween™ 20) according to the manufacturer's instructions. Samples were prepared in system buffer supplemented with 1 mg/mL of carboxymethyl-dextran (CMD, Sigma-Aldrich, code no. 86524). Samples were analyzed at 37°C. Approximately 12 000 RU goat anti-rabbit Fc fragment-specific antibody (GARbFcγ, Jackson ImmunoResearch Laboratories Inc, code no. 111-005-046) was immobilized according to the manufacturer's instructions using EDC/NHS coupling in both flow cells and in all eight channels. The sensor surface was finally saturated with 1 M ethanolamine. The binding kinetics of the respective rabbit mAbs against an HbA1c peptide analyte A01, relative molecular mass [M_r] of 1200) were generated and evaluated. The experimental assay setup is shown in Figure 1. Antibodies were injected at 200 nM concentration for 1 min at 10 μL/min. After capture, the flow rate was increased to 80 μL/min. Analyte A01 was injected for 3 min in a concentration series of 0 (buffer control), 3, 10, 30 (in duplicate), as well as 90 and 270 nM. The analyte dissociation was monitored for 5 min. After each analyte injection, the antibody capture system was fully regenerated by a 15 s injection of HBS-ET buffer at 20 μL/min; a 20 s injection at 20 μL/min with 10 mM glycine buffer pH 2.0, and two injections for 1 min at 20 μL/min with 10 mM glycine pH 2.25. Kinetic signatures were evaluated using a 1:1 binding model with local R_{max}.


Fig 1. Experimental assay setup for kinetic evaluation of six rabbit mAbs (A01 to A06) with specificity for human HbA1c peptide analyte A01. GARbFcγ was immobilized in flow cells 1 and 2 and in all eight channels. The rabbit mAbs were captured on GARbFcγ in flow cell 2. Human HbA1c was injected in both cells. Rabbit mAb X with unrelated target binding function and rabbit normal IgG1G2 were used as controls.

Binding specificity of rabbit mAb to various human HbA1c peptide derivatives

After the kinetic experiments, mAb A was selected to determine specificity against seven human HbA1c peptide analytes with M_r between 1100 and 1500 and a buffer control (Fig 2). The peptide analytes differed in their amino acid composition, sequence length, glycosylation, or N-/C-terminal chemical modifications. For example, A01 is glycosylated HbA1c positive control peptide and A02 and A05 are nonglycosylated HbA1c peptide negative controls. A07 is the HbA1A2 cross-reactive negative control. A03, A04, and A06 are positive controls for HbA1c peptide derivatives with chemical modifications or prolonged peptide sequences to avoid N- or C-terminal antibody peptide binding. Analytes A01, A03, A04, A05, and A06 were injected in a concentration series as previously described. Analytes A02 and A07 were injected at a higher concentration series from 0, 10, 30, 90 nM (in duplicate), and 270 nM to 810 nM to better identify unwanted cross-reactive binding.


Fig 2. Experimental assay setup for the specificity test of mAb A against seven human HbA1c peptide analytes, A01 to A07. GARbFcγ was immobilized in flow cells 1 and 2 and in all eight channels. Rabbit mAb A was captured on GARbFcγ in flow cell 2. Human HbA1c peptide analyte was injected in both cells.

Kinetic Screening in the Antibody Development Process

Antibody Methods and Protocols pp 171-18

Chapter 11

Kinetic Screening in the Antibody Development Process

Michael Schröml and Matthias Biehl

Abstract

Kinetic screening is of paramount importance when it is to select custom-made antibodies, tailored for their respective scientific, diagnostic, or pharmaceutical application. Here a kinetic screening protocol is presented. The assay is based on an ELISA method using a single injection. The method enables the screening of antibodies with different kinetic signatures.

if life ($t_{1/2, \text{off}}$), Binding Late,

4) means to select antigenic kinetic rate properties, dissociation constant K_d (M), in a workflow, where it ELISA-based screening, information about an antigen is of potentially suitable

asson resonance (SPR) uses six easy-to-access pictures to quickly select sample sets. The screening captured on the sensor binding signal at the end

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11 Kinetic Screening in the Antibody Development Process 177


Fig. 3. Exemplary valence analysis of the data set from Fig. 2. Binding Late (RU) is plotted over the antibody Capture Level (RU). Corridors are formed by lines, indicating areas of the same Molar Ratios. The lines help to quickly classify antibodies according to their valences. The antibodies populate all corridors and therefore bind to their 72-kDa antigen with different valences. Black triangles indicate the selected clone cultures from Fig. 2.

To graphically visualize the valence analysis, Binding Late is plotted over the antibody Capture Level. Valence corridors, calculated using the MR formula, allocate the antibodies according to their virtual Binding Late values (Fig. 3), e.g., for the plotting of the MR=2 trendline use the formula: BL at MR(2) = (MW(antigen) \times 2 \times Capture Level (RU)) / MW (antibody)). Just replace the Molar Ratio values to calculate different trendlines.

When the antibody shows a MR=1, it binds to an antigen with single valence. Obviously, there is some steric hindrance, which avoids full bivalent antibody binding.

When the antibody shows MR=2, it is able to simultaneously bind to two antigens.

Here antibodies in the range of MR=1 and MR=2 were selected, when they showed a suitable BL/SL ratio (Fig. 3) at the same time. For further details about the valence analysis see Notes 9 and 11.


3.5. Data Analysis for k_d Ranking

If possible, calculate the "off rate", which means the antigen dissociation rate constant k_d (1/s) (6). The dissociation rate constant k_d describes the decay of complexes per second. k_d is calculated from a linear regression fitting of the antigen dissociation phase, which is a function of the Biaevaluation software. Using the dissociation rate constant, the complex half life $t_{1/2, \text{off}}$ in minutes

Epitope binning analysis for further differentiation of therapeutic antibodies and diagnostic reagents

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Epitope binning analysis for further differentiation of therapeutic antibodies and diagnostic reagents
 Linnea Nygren-Baboi, Michael Murphy, and Anna Jansson
 GE Healthcare Bio-Sciences AB, SE-751 84 Uppsala, Sweden

Background

Our interest with this study was to establish a well characterized and easily accessible model system to be used for epitope binning assays using Biacore® and/or surface plasmon resonance (SPR) systems. Our criteria for the model system were that a monoclonal antigen and at least two antibodies should bind to the antigen simultaneously. The system we analyzed was epidermal growth factor receptor (EGFR) and three therapeutic antibodies targeting the extracellular domain of EGFR. The antibodies were panitumumab (Vectibix®), cetuximab (Erbitux®), and matuzumab (Cytuxen®). According to the previously published X-ray crystal structure, the Pan- fragment of matuzumab (Pan2000) interacts with an epitope adjacent to the cetuximab epitope and should thus be able to form a sandwich complex with cetuximab. Panitumumab and cetuximab, on the other hand, recognize the same epitope and should thus block each other. The major aim of the study was to evaluate the results from epitope binning using sandwich and premix assays and our interpretation of the findings.

Epitope binning EGFR (Sandwich)

Experimental conditions

- Antigen: epidermal growth factor (EGF), recombinant (EGF-EGF), and recombinant (EGF-EGF) in water (0.1M NaCl, 0.05M HEPES, pH 7.4)
- Antibody 1: Pan-EGFR
- Antibody 2: Cetuximab
- Antibody 3: Panitumumab
- Antibody concentrations: 4 nM, 40 nM, and 400 nM
- Flow rate: 30 µL/min
- Regeneration: 2 M NaOH, 10 min in buffer pH 2.0 (PharMA)

Results

The binding was observed in all of the combinations tested. The results from the different assays are shown in Figure 1. The results show that Pan-EGFR and Cetuximab can form a sandwich complex, while Panitumumab and Cetuximab cannot. This was not observed for Panitumumab.

Epidermal growth factor receptor (EGFR)

EGFR is a transmembrane protein and a member of the ErbB family of receptors, a subfamily of other closely related receptor tyrosine kinases: EGFR (ErbB-1, HER2/neu (ErbB-2), Her 3 (ErbB-3), and Her 4 (ErbB-4). Mutations affecting EGFR expression or activity could result in different cancer types. Anti-EGFR monoclonal antibodies for therapeutic use block the extracellular domain of EGFR and block ligand-induced EGFR tyrosine kinase activation. In this study, we used the extracellular domain (EGF-EGF) of human EGFR (ErbB1 domain).

Epitope binning EGFR (Premix)

Experimental conditions

- Antigen: recombinant EGF (EGF-EGF), recombinant (EGF-EGF), recombinant (EGF-EGF), recombinant (EGF-EGF) in water (0.1M NaCl, 0.05M HEPES, pH 7.4)
- Antibody 1: Pan-EGFR
- Antibody 2: Cetuximab
- Antibody 3: Panitumumab
- Antibody concentrations: 4 nM, 40 nM, and 400 nM
- Flow rate: 30 µL/min
- Regeneration: 2 M NaOH, 10 min in buffer pH 2.0 (PharMA)

Results

The results show that all three of the antibodies were blocking each other.

Sandwich assay

In this assay, the first antibody is immobilized on the sensor chip. The second antibody is then injected and binds to the antigen. The third antibody is then injected and binds to the antigen. The results show that Pan-EGFR and Cetuximab can form a sandwich complex, while Panitumumab and Cetuximab cannot.

Premix assay

In this assay, the antigen is immobilized on the sensor chip. The antibodies are then injected and bind to the antigen. The results show that all three of the antibodies were blocking each other.

Fab-EGFR complex

X-ray crystallographic data of the Fab fragments of cetuximab (C190) and matuzumab (C200) in complex with the EGF extracellular domain (EGF-EGF) are shown. The epitopes are shown in red and blue. The results show that Pan-EGFR and Cetuximab can form a sandwich complex, while Panitumumab and Cetuximab cannot.

Results presented as heat map and bin chart

The heat map shows the results from the different assays. The bin chart shows the results from the different assays. The results show that Pan-EGFR and Cetuximab can form a sandwich complex, while Panitumumab and Cetuximab cannot.

References

- Schmidt, L. et al. Molecular binding of EGFR presents the conformational requirements for the extracellular domain of EGFR. *PLoS ONE* 10(12): e0183888 (2017).
- Li, L. et al. Structural basis for inhibition of epidermal growth factor receptor by panitumumab. *Chem Biol* 16(11): 1201-1211 (2013).
- Schmidt, L. et al. Epitope binning of antibodies: conformational requirements for the extracellular domain of EGFR. *PLoS ONE* 10(12): e0183888 (2017).

Conclusions

Combination of data from different X-ray structures may not represent the "real life" scenario. The results from the different assays show that Pan-EGFR and Cetuximab can form a sandwich complex, while Panitumumab and Cetuximab cannot. This was not observed for Panitumumab. The model system did not fulfill our criteria as all antibodies recognized the same epitope or shared a common epitope. Ability to measure binding at low immobilization levels and low sample concentrations requires a sensitive instrument — in this case, Biacore S—no distinguishable binding from unbound antibodies. The response may occur when working with high surface and sample concentrations.

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