


Cytiva Webinar

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

- haemagglutinin (HA)
- Y anti-influenza serum

Biacore™ biosensor assays for quantitation of influenza virus and HCP

Biacore biosensor technology is used for virus haemagglutinin (HA) and host cell proteins (HCPs) in influenza vaccines. Accurate quantitation of influenza virus is achieved using an inhibition assay. The HA method shows high precision, and recovery as compared to immunodiffusion (SRID) which is the mainstay today. In addition, the analysis time is significantly shorter than the SRID method. Assays performed on Biacore systems demonstrate higher sensitivity and specificity compared to the SRID method, as a cost-effective replacement of existing methodologies.

Introduction

In vaccine production, the quantitation of concentration and HCP impurities is critical for formulation. Most influenza vaccines are produced in fertilized eggs, but the fear of a new pandemic has also increased the demand for short cycle systems. The need for shorter cycles has also increased the demand for analytical tools. In this Application note, the quantitation of influenza virus (I) and HCP technology are described. The potency of influenza vaccines is determined by quantitation of HA using indirect methods. The most common method used is SRID where vaccine dilution is diffused through a gel containing a specific antigen. SRID is a labor-intensive method and has disadvantages including low precision and

Biacore influenza virus quantitation assay

Recombinant HA antigen was immobilized on the dextran matrix to a level of 4 000 to 10 000 response units (RU) with standard amine coupling (Fig 1). Sera-containing antibodies with affinity for the specific influenza strains were diluted to give a response of about 1000 RU for a 400 s injection over the surface. This dilution factor resulted in a dynamic range of the standard curves of 0.5 to 8 µg virus/ml (Fig 1C). To stabilize the assay, 5 to 10 start-up cycles were run first with serum only, followed by regeneration.



Fig 1. Inhibition assay principle. (A) HA is first immobilized on the dextran matrix (red circles). Virus is then mixed with a fixed concentration of serum and injected over the surface. Free antibodies (not bound to virus at equilibrium) bind to the surface HA, giving a response. Low concentration of virus in the sample (A) gives high antibody binding, while high virus concentration (B) results in low binding level. (C) Overlay plot showing sensorgrams of injected serum mixed with a concentration series of virus standard. Report points (marked as X) were taken before and after injection to measure response levels. The surface was regenerated after each injection to prepare for a new sample.

Samples were analyzed without pretreatment. Mixtures of virus reference antigen or a virus sample with unknown concentration and a fixed dilution of anti-influenza serum were injected for 400 s (Fig 1A and 1B), followed by regeneration of the surface using freshly prepared 50 mM HCl, 0.05% Surfactant P20 for 30 s. Process samples were diluted so that at least two dilutions per sample reached a concentration within the standard curve range. Standard curves were run first, in the middle, and last in the assays. An interpolated calibration, performed by the software, was used to ensure high precision in calculated concentrations.

In the trivalent assay, three recombinant HA proteins (A/H1N1, A/H3N2, and B) were immobilized in three different flow cells on the same sensor chip. The three corresponding reference antigens and sera were mixed to form a trivalent standard solution. The standards and samples were allowed to pass over the surfaces, analyzing all three strains simultaneously.

Biacore biosensor HCP assay

HCP from MDCK cells were quantitated using a direct binding assay. Affinity-purified rabbit antibodies against MDCK cell lysate were immobilized using standard amine coupling to a response of 7 000 RU. The samples were diluted in HBS-EP+ with 1% BSA and injected over the surface for 300 s. The surface was regenerated with 50 mM NaOH. Standard curves using MDCK cell lysate were run first, in the middle, and last in measurement series. An interpolated calibration, performed by the software, was used to ensure high precision in calculated concentrations.

SRID assay, human influenza virus

A small portion of each serum was labeled with Cy3. Unlabeled serum and Cy3-labeled serum (< 0.5% of the total serum amount; did not interfere with the precipitation) were mixed with molten agarose (1%) and cast in a mould creating holes in the gel. Samples and reference antigens were treated with Zwittergent (1%) for 30 min at -22°C and added to the holes in the gel. The gels were incubated for 15 to 18 h at -22°C, and subsequently dried and scanned at 450 nm in Typhoon. Evaluation was performed by comparing ring areas of samples with reference antigens using ImageQuant TL software.

Bradford total protein assay

The method was performed according to the manufacturer's instructions. Bradford solution (200 µl) was mixed with BSA standard or diluted samples (20 µl) in a 96-well plate format. The plate was incubated for ~1 h, and measured in a plate reader at 595 nm. The evaluation was performed in the Softmax™ Pro v5 software.


* This assay was developed as part of an MSc thesis by E. Anderson, Uppsala University, 2009.

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
Karlsson, R. Biophysical Reviews 1-12 (2016)
Pol et.al. Analytical Biochemistry 510, 88-97 (2016)

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
Assessment of protein concentration using a calibration-free method

Kevin Lindquist, Åsa Sporang, Robert Karlsson, Anna Sylvan, Anita Larsson, and Ewa Pol. GE Healthcare, Uppsala, Sweden.


Introduction

In development and production of biotherapeutics it is important to determine the concentration of functionally active molecules in a given sample. Many established methods for protein concentration measurement do not distinguish between active and inactive molecules, or rely on comparison with a standard that may or may not be a mimic of the actual sample. The method discussed here, Calibration-Free Concentration Analysis (CFCA), provides a highly sensitive and robust biosensor-based determination of the protein concentration and does not require a standard. The method relies on measurement of analyte binding to a target immobilized on a sensor surface at varying flow rates, under conditions where the observed rate of binding is partially or completely limited by transport of analyte molecules to the sensor surface. This transport is a diffusion controlled process that depends on the analyte concentration. The concentration is obtained by running Biacore® binding experiments at different flow rates and fitting the binding data to the model describing the process.

In Biacore instruments samples are injected in a micro flow system and transported in laminar flow to the sensor surface. Molecules reach the sensor surface from bulk solution by a diffusion-controlled transport process. In addition to the concentration of analyte molecules, factors influencing the transport rate include the diffusion coefficient, flow cell dimensions, and flow rate. The balance between the transport rate and the binding rate determines whether the observed binding will be transport limited or reaction limited. For successful CFCA, the observed binding rate must be at least partially limited by transport. In this poster, we discuss how suitable flow rates, immobilization levels, and binding data are selected for successful implementation of CFCA.

How to recognize transport limited data?

- Reaction limited data
 - Binding rates are not affected by flow rate variation
 - Data not suitable for CFCA
- Transport limited data
 - Binding rates increase with increasing flow rate
 - Data suitable for CFCA
- Partially transport limited data
 - Initial binding is affected by flow rate variation
 - Late binding is less affected by the flow rate
 - Data potentially suitable for CFCA

Immobilization level and concentration range

Calibration-free concentration analysis of 2-fold dilution series of β_2 -microglobulin was performed on four surfaces immobilized to increasing levels of anti- β_2 -microglobulin antibody. Immobilization levels are expressed as response units, RU, divided by the molecular weight of interaction partner, kDa, to compensate for the difference in response between differently sized molecules. The protein solutions were injected at two flow rates, 5 and 90 μ l/min.

The concentration determined by the method is shown in each panel. Sensorgrams framed in green are obtained under experimental conditions that return the same calculated concentration after adjustment for the dilution factor. The concentration given by sensorgrams framed in red are overestimated.

The most apparent visual indicator of data suitable for CFCA is a significant separation of the binding curves obtained at high and low flow rates.



Flow rates

Concentrations were determined for two independent antibody model systems using binding data obtained at flow rates of 5, 10, 20, 40, 60, 80 and 100 μ l/min. The binding data was analyzed in all 21 pairwise combinations of flow rates, as well as in the whole data set of 7 flow rates. Results were normalized to the value obtained from 7 flow rates (shown as value 1 with thinner lines to indicate +/- SE) to allow comparison between the model systems.


The results indicate that reliable concentration data is obtained when using two flow rates that differ by a factor of 10 or more.

References

Glosser, K.W., 1993. Anal Biochem, 213: 152-161.
 Karlsson, R., Fägerlin, L., Malmberg, M., Persson, B., 1993. J. Immunol Methods, 168: 75-86.
 Karlsson, R., Bock, M., Fägerlin, L., Persson, B., 1994. Comparison Methods Diagram, 6: 99-110.
 Christman, L. H., 1997. Anal Biochem, 249: 313-316.
 Sjöqvist, K., Malmberg, M., Bock, M., Fägerlin, L., Persson, B., 2002. Biochemistry, 41: 8263-8276.
 Hyatt, G., Horton, M., Doyle, L., and Chabon, J.H., 1997. Biophys Chem, 64: 127-137.

Concentration of 8 commercial antibodies

The CFCA of eight commercial antibodies, denoted as A to H, was performed under transport limited conditions using anti-Fc antibody and protein A. The concentration values were compared to the concentration declared by manufacturers. Activity is given as the ratio of measured and specified concentrations, respectively. The concentrations of all antibodies were in good agreement when measured against two different interaction partners although could differ as much as 50% when compared to the concentration specified by manufacturers.

Conclusions

- Calibration-free concentration analysis on Biacore systems open new possibilities for concentration measurements that relate directly to determination of the concentration of a particular binding site. Measurements are direct, label-free, and do not require a calibration standard. Robust and reliable results are obtained when experiments are performed using:
 - at least two flow rates, where the ratio between flow rates is > 10
 - a high surface concentration, typically \geq 25 RU/kDa (immobilized level/ M_w of interaction partner)
 - protein dilutions that give concentrations \leq 30 nM
- In the example, the concentrations of eight antibodies were in good agreement when measured against two interaction partners targeting the antibody Fc part. However, concentrations determined by CFCA differed by up to 50% from the values specified, depending on the method used. More work is required to explain the discrepancies.

GE, the logo of GE and GE Healthcare are trademarks of General Electric Company. Biacore is a trademark of GE Healthcare Company. © 2008 GE Healthcare Company. All rights reserved.

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
【Tips】拡散係数の求めかた Diffusion Coefficient Calculator / Converter

You are here: [Biacore Life Sciences](#) > [Application Support](#) > [Laboratory Guidelines etc.](#) > [Diffusion Coefficient Calculator](#)

Training Portal

- Application Support
 - Laboratory Guidelines etc.
 - Training
 - Product Support & Service
- Web Services
 - My Account

← BACK PRINT

Diffusion Coefficient Calculator / Converter

This on-line tool is designed to help you calculate diffusion coefficients for use in Calibration-Free Concentration Analysis assays. It is accessible only via a valid product key associated with the appropriate types of Biacore software.

Note that the diffusion coefficient at 20 C is entered in the run method as a variable for CFCA determination. The software converts the value to the run temperature if necessary.

Calculate diffusion coefficient at 20°C

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