

## 仕様

Column	Tube inner diameter (mm)	Cross sectional area (cm <sup>2</sup> )	Tube height (mm)	Bed height (mm)		Bed vol. (ml)		Design pressure (MPa) *
				min.	max.	min.	max.	
FineLINE Pilot 35	35	10	330	30	150	29	144	2
FineLINE 70	70	38	350	30	150	115	577	2
FineLINE 70L	70	38	700	50	300	192	1,155	2
FineLINE 100P	100	79	350	30	150	236	1,178	2
FineLINE 100PL	100	79	700	50	300	393	2,355	2
FineLINE 200P	200	314	350	30	150	942	4,710	2
FineLINE 200PL	200	314	700	50	300	1,570	9,420	2
FineLINE 350P	350	962	350	30	150	2,884	14,424	2
FineLINE 350PL	350	962	700	50	300	4,808	28,848	2

\* 設計圧力は充填圧と操作圧の合計です。