

# HyClone classical media

CELL CULTURE

## Dulbecco's Modified Eagle's Medium - High Glucose Liquid media

Component description	Code number SH30019		Code number SH30022, SH30308		Code number SH30081, SH30315		Code number SH30243, SH30317	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
<b>Inorganic salts</b>								
Calcium chloride	200	1.8021	200	1.8021	200	1.8021	200	1.8021
Ferric nitrate-9H <sub>2</sub> O	0.1	0.0002	0.1	0.0002	0.1	0.0002	0.1	0.0002
Potassium chloride	400	5.3655	400	5.3655	400	5.3655	400	5.3655
Magnesium sulfate	97.67	0.8112	97.67	0.8112	97.67	0.8112	97.67	0.8112
Sodium chloride	6400	109.514	6400	109.514	6400	109.514	6400	109.514
Sodium phosphate monobasic H <sub>2</sub> O	125	0.9059	125	0.9059	125	0.9059	125	0.9059
<b>Amino acids</b>								
L-Arginine-HCl	84	0.3987	84	0.3987	84	0.3987	84	0.3987
L-Cystine-2HCl	62.57	0.1998	62.57	0.1998	62.57	0.1998	62.57	0.1998
L-Glutamine	584	3.9959	584	3.9959	0	0	584	3.9959
Glycine	30	0.3996	30	0.3996	30	0.3996	30	0.3996
L-Histidine-HCl-H <sub>2</sub> O	42	0.2004	42	0.2004	42	0.2004	42	0.2004
L-Isoleucine	104.8	0.7989	104.8	0.7989	104.8	0.7989	104.8	0.7989
L-Leucine	104.8	0.7990	104.8	0.7990	104.8	0.7990	104.8	0.7990
L-Lysine-HCl	146.2	0.8004	146.2	0.8004	146.2	0.8004	146.2	0.8004
L-Methionine	30	0.2011	30	0.2011	30	0.2011	30	0.2011
L-Phenylalanine	66	0.3995	66	0.3995	66	0.3995	66	0.3995
L-Serine	42	0.3997	42	0.3997	42	0.3997	42	0.3997
L-Threonine	95.2	0.7992	95.2	0.7992	95.2	0.7992	95.2	0.7992
L-Tryptophan	16	0.0783	16	0.0783	16	0.0783	16	0.0783
L-Tyrosine-2Na-2H <sub>2</sub> O	103.79	0.3974	103.79	0.3974	103.79	0.3974	103.79	0.3974
L-Valine	93.6	0.7990	93.6	0.7990	93.6	0.7990	93.6	0.7990

Component description	Code number SH30019		Code number SH30022, SH30308		Code number SH30081, SH30315		Code number SH30243, SH30317	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
<b>Vitamins</b>								
Calcium D-pantothenate	4	0.0084	4	0.0084	4	0.0084	4	0.0084
D-Pantothenic acid Na salt	0	0	0	0	0	0	0	0
Choline chloride	4	0.0286	4	0.0286	4	0.0286	4	0.0286
Folic acid	4	0.0091	4	0.0091	4	0.0091	4	0.0091
Myo-inositol	7	0.0389	7	0.0389	7	0.0389	7	0.0389
Niacinamide	4	0.0328	4	0.0328	4	0.0328	4	0.0328
Pyridoxine-HCl	4	0.0195	4	0.0195	4	0.0195	4	0.0195
Riboflavin	0.4	0.0011	0.4	0.0011	0.4	0.0011	0.4	0.0011
Thiamine-HCl	4	0.0119	4	0.0119	4	0.0119	4	0.0119
<b>Other</b>								
D-Glucose	4500	24.9778	4500	24.9778	4500	24.9778	4500	24.9778
HEPES	0	0	0	0	0	0	0	0
Phenol red-Na	15.9	0.0422	15.9	0.0422	15.9	0.0422	15.9	0.0422
Sodium pyruvate	0	0	0	0	0	0	110	0.9996
Sodium bicarbonate	3700	44.0424	3700	44.0424	3700	44.0424	3700	44.0424

## Dulbecco's Modified Eagle's Medium - High Glucose Liquid media

Component description	Code number SH30249		Code number SH30262, SH30319		Code number SH30284		Code number SH30285, SH30320	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
<b>Inorganic salts</b>								
Calcium chloride	200	1.8021	0	0	200	1.8021	200	1.8021
Ferric nitrate-9H <sub>2</sub> O	0.1	0.0002	0.1	0.0002	0.1	0.0002	0.1	0.0002
Potassium chloride	400	5.3655	400	5.3655	400	5.3655	400	5.3655
Magnesium sulfate	97.67	0.8112	0	0	97.67	0.8112	97.67	0.8112
Sodium chloride	4750	81.2799	6400	109.514	6400	109.514	6400	109.514
Sodium phosphate monobasic H <sub>2</sub> O	125	0.9059	125	0.9059	125	0.9059	125	0.9059
<b>Amino acids</b>								
L-Arginine-HCl	84	0.3987	84	0.3987	84	0.3987	84	0.3987
L-Cystine-2HCl	62.57	0.1998	62.57	0.1998	62.57	0.1998	62.57	0.1998
L-Glutamine	584	3.9959	0	0	584	3.9959	0	0
Glycine	30	0.3996	30	0.3996	30	0.3996	30	0.3996
L-Histidine-HCl-H <sub>2</sub> O	42	0.2004	42	0.2004	42	0.2004	42	0.2004
L-Isoleucine	104.8	0.7989	104.8	0.7989	104.8	0.7989	104.8	0.7989
L-Leucine	104.8	0.7990	104.8	0.7990	104.8	0.7990	104.8	0.7990
L-Lysine-HCl	146.2	0.8004	146.2	0.8004	146.2	0.8004	146.2	0.8004
L-Methionine	30	0.2011	30	0.2011	30	0.2011	30	0.2011
L-Phenylalanine	66	0.3995	66	0.3995	66	0.3995	66	0.3995
L-Serine	42	0.3997	42	0.3997	42	0.3997	42	0.3997
L-Threonine	95.2	0.7992	95.2	0.7992	95.2	0.7992	95.2	0.7992
L-Tryptophan	16	0.0783	16	0.0783	16	0.0783	16	0.0783
L-Tyrosine-2Na-2H <sub>2</sub> O	103.79	0.3974	103.79	0.3974	103.79	0.3974	103.79	0.3974
L-Valine	93.6	0.7990	93.6	0.7990	93.6	0.7990	93.6	0.7990
<b>Vitamins</b>								
Calcium D-pantothenate	4	0.0084	0	0	4	0.0084	4	0.0084
D-Pantothenic acid Na salt	0	0	4	0.0166	0	0	0	0
Choline chloride	4	0.0286	4	0.0286	4	0.0286	4	0.0286
Folic acid	4	0.0091	4	0.0091	4	0.0091	4	0.0091
Myo-inositol	7	0.0389	7	0.0389	7	0.0389	7	0.0389
Niacinamide	4	0.0328	4	0.0328	4	0.0328	4	0.0328
Pyridoxine-HCl	4	0.0195	4	0.0195	4	0.0195	4	0.0195
Riboflavin	0.4	0.0011	0.4	0.0011	0.4	0.0011	0.4	0.0011
Thiamine-HCl	4	0.0119	4	0.0119	4	0.0119	4	0.0119
<b>Other</b>								
D-Glucose	4500	24.9778	4500	24.9778	4500	24.9778	4500	24.9778
HEPES	5958	25.002	0	0	0	0	0	0
Phenol red-Na	15.9	0.0422	15.9	0.0422	0	0	15.9	0.0422
Sodium pyruvate	0	0	0	0	0	0	110	0.9996
Sodium bicarbonate	3700	44.0424	3700	44.0424	3700	44.0424	3700	44.0424

## Dulbecco's Modified Eagle's Medium - High Glucose Liquid media

Component description	Code number SH30563		Code number SH30585		Code number SH30604		Code number SH30606	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
<b>Inorganic salts</b>								
Calcium chloride	200	1.8021	200	1.8021	200	1.8021	200	1.8021
Ferric nitrate-9H <sub>2</sub> O	0.1	0.0002	0.1	0.0002	0.1	0.0002	0.1	0.0002
Potassium chloride	400	5.3655	400	5.3655	400	5.3655	400	5.3655
Magnesium sulfate	97.67	0.8112	97.67	0.8112	97.67	0.8112	97.67	0.8112
Sodium chloride	6400	109.514	6400	109.514	6400	109.514	6400	109.514
Sodium phosphate monobasic H <sub>2</sub> O	125	0.9059	125	0.9059	125	0.9059	125	0.9059
<b>Amino acids</b>								
L-Arginine-HCl	84	0.3987	84	0.3987	84	0.3987	84	0.3987
L-Cystine-2HCl	62.57	0.1998	62.57	0.1998	62.57	0.1998	0	0
L-Glutamine	0	0	0	0	0	0	0	0
Glycine	30	0.3996	30	0.3996	30	0.3996	30	0.3996
L-Histidine-HCl-H <sub>2</sub> O	42	0.2004	42	0.2004	42	0.2004	42	0.2004
L-Isoleucine	104.8	0.7989	104.8	0.7989	104.8	0.7989	104.8	0.7989
L-Leucine	104.8	0.7990	104.8	0.7990	104.8	0.7990	104.8	0.7990
L-Lysine-HCl	146.2	0.8004	146.2	0.8004	146.2	0.8004	146.2	0.8004
L-Methionine	30	0.2011	30	0.2011	30	0.2011	0	0
L-Phenylalanine	66	0.3995	66	0.3995	66	0.3995	66	0.3995
L-Serine	42	0.3997	42	0.3997	42	0.3997	42	0.3997
L-Threonine	95.2	0.7992	95.2	0.7992	95.2	0.7992	95.2	0.7992
L-Tryptophan	16	0.0783	16	0.0783	16	0.0783	16	0.0783
L-Tyrosine-2Na-2H <sub>2</sub> O	103.79	0.3974	103.79	0.3974	103.79	0.3974	103.79	0.3974
L-Valine	93.6	0.7990	93.6	0.7990	93.6	0.7990	93.6	0.7990
<b>Vitamins</b>								
Calcium D-pantothenate	4	0.0084	4	0.0084	4	0.0084	4	0.0084
D-Pantothenic acid Na salt	0	0	0	0	0	0	0	0
Choline chloride	4	0.0286	4	0.0286	4	0.0286	4	0.0286
Folic acid	4	0.0091	4	0.0091	4	0.0091	4	0.0091
Myo-inositol	7	0.0389	7	0.0389	7	0.0389	7	0.0389
Niacinamide	4	0.0328	4	0.0328	4	0.0328	4	0.0328
Pyridoxine-HCl	4	0.0195	4	0.0195	4	0.0195	4	0.0195
Riboflavin	0.4	0.0011	0.4	0.0011	0.4	0.0011	0.4	0.0011
Thiamine-HCl	4	0.0119	4	0.0119	4	0.0119	4	0.0119
<b>Other</b>								
D-Glucose	4500	24.9778	4500	24.9778	4500	24.9778	4500	24.9778
HEPES	5958	25.002	0	0	0	0	0	0
Phenol red-Na	15.9	0.0422	0	0	0	0	15.9	0.0422
Sodium pyruvate	110	0.9996	0	0	110	0.9996	110	0.9996
Sodium bicarbonate	3700	44.0424	3700	44.0424	3700	44.0424	3700	44.0424

## Dulbecco's Modified Eagle's Medium - High Glucose Liquid media

Code number  
SH30607

Component description	mg/L	mmol/L
<b>Inorganic salts</b>		
Calcium chloride	200	1.8021
Ferric nitrate-9H <sub>2</sub> O	0.1	0.0002
Potassium chloride	400	5.3655
Magnesium sulfate	97.67	0.8112
Sodium chloride	6400	109.514
Sodium phosphate monobasic H <sub>2</sub> O	0	0
<b>Amino acids</b>		
L-Arginine-HCl	84	0.3987
L-Cystine-2HCl	62.57	0.1998
L-Glutamine	0	0
Glycine	30	0.3996
L-Histidine-HCl-H <sub>2</sub> O	42	0.2004
L-Isoleucine	104.8	0.7989
L-Leucine	104.8	0.7990
L-Lysine-HCl	146.2	0.8004
L-Methionine	30	0.2011
L-Phenylalanine	66	0.3995
L-Serine	42	0.3997
L-Threonine	95.2	0.7992
L-Tryptophan	16	0.0783
L-Tyrosine-2Na-2H <sub>2</sub> O	103.79	0.3974
L-Valine	93.6	0.7990
<b>Vitamins</b>		
Calcium D-pantothenate	4	0.0084
D-Pantothenic acid Na salt	0	0
Choline chloride	4	0.0286
Folic acid	4	0.0091
Myo-inositol	7	0.0389
Niacinamide	4	0.0328
Pyridoxine-HCl	4	0.0195
Riboflavin	0.4	0.0011
Thiamine-HCl	4	0.0119
<b>Other</b>		
D-Glucose	4500	24.9778
HEPES	0	0
Phenol red-Na	15.9	0.0422
Sodium pyruvate	110	0.9996
Sodium bicarbonate	3700	44.0424

### [cytiva.com/hyclone](https://www.cytiva.com/hyclone)

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. HyClone is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY13878-10Aug20-DF

