



## HiTrap™ NHS-activated HP column

Dear customer,

Dryness can sometimes be observed in the prepacked columns. Please note that the dryness does not affect the performance of the columns and will disappear during column preparation.

[cytiva.com/hitrap](https://www.cytiva.com/hitrap)

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva.

HiTrap is trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2022 Cytiva

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY26994-15Feb22-FL